Superior

Anti-Leak Face Mask

& accessories against PM2.5



MOTEX	face mask and	l accessory	facial fit	comparison
-------	---------------	-------------	------------	------------

Diamond Shape Face Mask-Combination	*Average Facial Fit	Average upgrade	Diamond Shape Face Mask-Combination	*Average Facial Fit	Average upgrade
Anti-Leak Diamond Shape Face Mask + Loop-Adjuster +Leak-proof sheet	30	800%	Anti-Leak Flat Shape Face Mask + Loop-Adjuster +Leak-proof sheet	20	500%
Anti-Leak Diamond Shape Face Mask + Loop-Adjuster	22	600%	Anti-Leak Flat Shape Face Mask + Loop-Adjuster	16	400%
Anti-Leak Diamond Shape Face Mask	16	400%	Anti-Leak Flat Shape Face Mask	10	200%
Diamond Shape Face Mask	10	200%	Flat Shape Face Mask	%3.1	

* Quantitative Fit Test Fit Factor C_i : the concentration measured outside the respirator C_0 : the concentration measured inside the respirator

According to the protective masks (N95/FFP2) D2 facial fit testing methods in CNS14774 and CNS14755, the facial fit of face mask must reach 100. According to the 100-person test of the Institute of Occupational Safety and Health of Taiwan (Taiwan IOSH), the facial fit of random sampling two types of N95 masks in the market is 13.9 and 21.8 respectively, and the facial fit of traditional flat shape surgical masks is only 3.1. (Data extracted from the Taiwan IOSH website)



Motex Healthcare Corp.

19F-3, No.508, Sec. 5, Chung Hsiao East Rd. Taipei 110, Taiwan Tel:886-2-2727-7347 Fax:886-2-2727-8323 E-mail:info@motex com http://www.motex.com

